

## MULTI-WAVELENGTH ELLIPSOMETERS

Innovative ■ Powerful ■ Affordable

Film Sense Multi-Wavelength Ellipsometers use long-life LED's and a no-moving-parts ellipsometric detector to provide fast and reliable thin film measurements in an easy-to-use, compact system.

The film thickness and index of refraction of most transparent thin films can be determined with excellent precision and accuracy by a simple 1 second measurement. Optical constants  $n$  &  $k$  and other film properties can also be measured for many samples.

Multi-Wavelength Ellipsometry provides powerful thin film measurement capabilities, while at the price point of single wavelength ellipsometer and spectroscopic reflectometer systems. Film Sense ellipsometers are ideal for measurements in the research lab, classroom, in situ process chambers, industrial quality control, and more.


# What is Ellipsometry?

Ellipsometry measures the change in polarization state for light reflected from a Sample. The ellipsometric measurement is quantified by the formula below, where  $\rho$  is the complex ratio of the reflectivities for p- and s- polarized light ( $R_p$  and  $R_s$ ). The ellipsometric  $\Psi$  parameter is related to the magnitude of the complex ratio, and the ellipsometric  $\Delta$  parameter is the phase of the complex ratio.

$$\rho = \frac{R_p}{R_s} = \tan(\Psi) \cdot e^{i \cdot \Delta}$$


## Advantages of Ellipsometry

- Ellipsometry measures a ratio, so it is not sensitive to changes in the beam intensity, or sample imperfections that scatter light away.
- Since ellipsometry measures 2 quantities ( $\Psi$  and  $\Delta$ ) at each wavelength, it can determine 2 quantities, such as film thickness & index, or substrate n & k.
- The ellipsometric  $\Delta$  parameter is extremely sensitive to thin films, enabling accurate film thickness measurements down to 0 Å.



## Advantages of Multiple Wavelengths

- Enable unique determination of film thickness for transparent films (no thickness periodicity issues).
- Determine additional sample parameters, such as: surface roughness, multiple film thicknesses, index dispersion.
- Provide a consistency check on the data analysis—a “good” analysis model should fit the data at all measured wavelengths.
- For very thin films (< 20 nm), a multi-wavelength ellipsometric data set can provide information content similar to spectroscopic ellipsometric data (contact Film Sense and request our “Multi-Wavelength vs. Spectroscopic Data” white paper for more details).



*To determine sample parameters of interest, such as film thickness and index of refraction, an optical model is used to analyze the ellipsometric data.*

## Film Sense Multi-Wavelength Ellipsometer Technology

### KEY FEATURES

Multiple LED sources (either 4 or 8, with wavelengths ranging from 370 – 950 nm, depending on the system)

No moving parts in the ellipsometric detector\*

Excellent thickness precision, better than 0.0004 nm for many samples (for a 1 second acquisition), even for sub-monolayer film thicknesses

Integrated computer for instrument control and data analysis, with a web browser interface accessible from any modern computer, laptop, or tablet

Completely self-contained system

### BENEFITS

Long lifetimes (>50,000 hours), with no costly lamp changes, time consuming alignments or PM procedures

Fast measurement times (ellipsometric data in 1.7 ms, with the new *Fast Mode* feature) and long term reliability

Measurement precision that is only possible with an ellipsometer

No complicated software setup and maintenance

No external electronics box or fiber connections

**NOW AVAILABLE**

# Film Sense Multi-Wavelength Ellipsometer Systems

The 4th generation Film Sense Multi-Wavelength Ellipsometer systems are now available! The primary improvements in this generation are the additional wavelengths and spectral range for the FS-8 model, which further enhance the measurement capabilities over a wide range of thin film applications.



## Gen. 4 Common Specifications

- Compact optics: Source 125 x 80 x 60 mm  
Detector 110 x 80 x 60 mm
- Simple connections: +12V wall plug power supply, Ethernet, and Source-Detector link cables
- Motorized Source Polarizer:
  - provides automated instrument calibration
  - enables zone-averaged measurements, for improved measurement accuracy
- 4x more intensity (compared to original FS-1), and updated detector electronics provide improved measurement precision: 2x for ex situ, 4x for in situ

## FS-4

- 4 wavelengths, 450 – 660 nm spectral range (replaces the original FS-1)
- Excellent choice for measuring single layer transparent films in the 0 – 2  $\mu\text{m}$  thickness range, with precision down to 0.0004 nm

## FS-8

- 8 wavelengths, 370 – 950 nm spectral range
- The UV wavelength (370 nm) provides enhanced sensitivity when measuring very thin films < 10 nm
- The 3 longer wavelengths (735, 850, and 950 nm) enable the measurement of thicker transparent films (up to 5  $\mu\text{m}$ ), and absorbing semiconductor films (such as poly-Si, SiGe, amorphous-Si, etc.)
- Film resistivity measurements (using the Drude model) are also improved with the 3 longer wavelengths.
- 8 wavelengths and wider spectral range provide enhanced measurement capability for multilayer film stacks

## Film Sense Wavelengths

System	Number of Wvls	Wavelength (nm)								
		370	405	450	525	595	660	735	850	950
FS-1 (gen. 2 & 3)	4			X	X	X	X			
FS-1EX (gen. 2 & 3)	6		X	X	X		X		X	X
FS-4 (gen. 4)	4			X	X	X	X			
FS-8 (gen. 4)	8	X		X	X	X	X	X	X	X

\*US Patent #9,354,118



# Capabilities and Performance

Film Sense Multi-Wavelength ellipsometers excel at measuring the thickness and index of refraction of transparent single films. The upper thickness limit depends on the ellipsometer system (typically 2 – 5  $\mu\text{m}$ ), but is also dependent on the substrate and film optical constants. As with any ellipsometer system, a minimum film thickness (typically 10 nm) is required to obtain accurate index of refraction measurements (the new Liquid Cell option can remove this limitation).

Optically absorbing films can also be measured, but the data analysis becomes more complicated as the film optical constants (both n and k values) are required. The Film Sense software contains multiple methods for determining n&k values: 1) multi-sample analysis, 2) combined ellipsometry + transmission measurements, 3) immersion ellipsometry with the new Liquid Cell option, and 4) dispersion models. The upper thickness limit for absorbing films is strongly dependent on the type of material; for metallic films, the upper limit is typically 50 nm.

Multi-Wavelength Ellipsometry can also be used to measure multilayer film stacks (in some cases up to 5 layers), depending on the thicknesses and indices of refraction of the layers. Simulations can be performed in the Film Sense software to determine if a particular sample structure is possible. For some samples, surface roughness and index gradients in the film can also be characterized.

The typical Film Sense ellipsometer measurement Accuracy and Precision for a variety of samples, including a multi-layer sample, is shown in the table below. For more details on the testing methodology, contact Film Sense and request our "FS-8 Performance" white paper.

Sample	Parameter	Accuracy	Precision
#1. 2nm Native Oxide on Si	Thickness	0.013 nm	0.00026 nm
#2. 90nm Oxide on Si	Thickness	0.11 nm	0.00041 nm
	Index @ 633 nm	0.0005	7.3E-06
#3. 1000nm Oxide on Si	Thickness	0.60 nm	0.0038 nm
	Index @ 633 nm	0.0001	3.3E-06
#4. 100-50-100 nm ONO on Si	Top SiO <sub>2</sub> Thickness	0.04 nm	0.0023
	Si <sub>3</sub> N <sub>4</sub> Thickness	0.08 nm	0.0032
	Bottom SiO <sub>2</sub> Thickness	0.10 nm	0.0038
	SiO <sub>2</sub> Index @ 633 nm	0.0003	4.3E-06
#5. 6 nm TiO <sub>2</sub> on Si	Thickness	0.002 nm	0.00030 nm
	Index @ 633 nm	0.004	3.3E-04
#6. 70 nm Al <sub>2</sub> O <sub>3</sub> on Si	Thickness	0.042 nm	0.0010 nm
	Index @ 633 nm	0.0011	1.1E-05
	k @ 633 nm	0.0006	8.3E-06
#7. 500 nm SiN on Si	Thickness	0.59 nm	0.0036 nm
	Surface Roughness	0.004 nm	0.0011 nm
	Index @ 633 nm	0.0013	1.3E-05
	k @ 633 nm	0.0001	1.6E-06
#8. 300 nm Ta <sub>2</sub> O <sub>5</sub> on Glass	Thickness	0.29 nm	0.0031 nm
	Surface Roughness	0.12 nm	0.0012 nm
	Index @ 633 nm	0.0006	1.7E-05
#9. 60 nm a-Si on Glass	Thickness	0.08 nm	0.0009 nm
	Oxide Thickness	0.02 nm	0.00054 nm
	Band Gap	0.0024 eV	5.7E-05
#10. 10 nm TiN on Si	Thickness	0.042 nm	0.00047 nm
	Resistivity	1.1 uOhm-cm	0.020 uOhm-cm

## SEND US YOUR SAMPLES!

As thin film applications are so varied and diverse, the best way to determine if a Film Sense Multi-Wavelength Ellipsometer is right for your application is to perform demonstration measurements on your actual samples. Please contact us to discuss your application, and arrange for sample measurements.

## Standard Ex Situ Configuration

- 65° Angle of Incidence.
- Manual sample loading and height adjustment.
- Sample sizes up to 200 mm dia. and 20 mm in thickness.
- Sample tilt with +/-2° range
- Beam size on sample: 4 x 9 mm.
- Compact footprint (180 x 400 mm) and light (5 kg).



## Focused Beam Option


- Reduces beam size on sample to 0.8 x 1.9 mm or 0.25 x 0.55 mm.





# In Situ Features

Film Sense Multi-Wavelength Ellipsometers are ideal for in situ realtime monitoring and control of thin film deposition and etching processes.

- Completely self-contained system, with no external electronics box or fiber connections required
- LED light sources and no moving parts detector, for robust and reliable operation, and fast measurements
- Compact and light weight source and detector units (< 1 kg each)
- Optional adapters for mounting to standard 2.75" or 1.33" conflat flanges, with easy to adjust coarse and fine tilt stages
- Powerful software features for visualizing and analyzing dynamic ellipsometric data, with the new **Model Builder** feature to automate the analysis



FS-1 Mounted on AJA Sputter Chamber




FS-8 with Auto-Align option, mounted on Kurt Lesker ALD Chamber

# In Situ Capabilities

- Sub-monolayer thickness precision
- Determine film optical constants n&k and deposition rates, at multiple process conditions, without breaking vacuum
- Monitor and control the deposition of multilayer film structures
- FS-API interface for external software control (LabVIEW™ compatible)
- Applicable to most thin film deposition and etching techniques: Sputtering, ALD, ALE, MBE, CVD, PLD, etc.
- **Auto-Align** option provides automated in situ alignment of the ellipsometer Src and Det units (see <https://film-sense.com/film-sense-in-situ-monitoring-system/#auto-align> for a video demo)
- **Fast Mode** feature enables ultra fast (1.7 ms) ellipsometric data acquisition, while still providing sub-Å thickness precision

NEW



Data courtesy of the Army Research Lab







# FS-RT300 Automated Mapping System

Combines an FS-8 Multi-Wavelength Ellipsometer with a compact automated mapping stage to provide fast, accurate, and reliable film thickness uniformity measurements across a wafer.

## Features and Specifications

- 8 wavelengths of ellipsometric data (370, 450, 525, 595, 660, 735, 850, 950 nm), with long life LED sources, and no moving parts detector
- Accurate thickness measurements for most transparent thin films from 0 – 5  $\mu\text{m}$
- Typical time for wafer map: 2.5 minutes (49 points on a 300 mm diameter wafer)
- Typical thickness repeatability: 0.002 nm
- Integrated focusing probes, standard spot size: 0.8 x 1.9 mm (other spot sizes available)
- Compact footprint: 400x500 mm, 22 kg
- Stage travel: R (linear) 150 mm, resolution: 12  $\mu\text{m}$ , Theta (rotation) 360°, resolution: 0.05°
- USB camera option available
- Motorized Z-stage for sample auto alignment
- Contour and 3D plots of measured parameters
- Flexible Scan Pattern Editor



# Software

The Film Sense software acquires and analyzes the ellipsometric data, and reports the sample parameters (thickness, index of refraction, etc.) that are derived from the measurement. The Film Sense software runs on a computer which is inside the Detector unit, and a standard web browser provides the user interface for the software. Any desktop, laptop, or tablet computer that supports a modern web browser (Windows, Mac OS X, Linux, iOS, Android) can operate the Film Sense ellipsometer using its Ethernet connection (no Internet or web access is required). A major advantage to the web browser interface is that no software installation is required, which greatly simplifies the setup and operation of Film Sense ellipsometers. (An Analysis-Only version of the software is also available to support offline data analysis.)

The **Single Measurement** screen makes routine sample measurements as easy as clicking a button.

Fit_Diff	0.0016
Thick(nm).1	70.584
n(633nm).1	1.6434
n_slope.1	0.0061911

The new **Model Builder** feature automates the data analysis process, and provides expert advice on the best model (see <https://film-sense.com/software/#Model-Builder> for a video demo)

Description	Fit_Diff	Thick Err.	Thick(nm).1
Book Value Opt. Const.	0.0506	0.9773	81.45
with Surface Roughness	0.0497	1.383	81.038
Ideal Cauchy Model	0.0303	2.39	87.089
Cauchy with Roughness	0.0282	2.627	86.094
<b>Cauchy with Roughness+Grade (BEST)</b>	<b>0.0024</b>	<b>0.2407</b>	<b>87.119</b>
Cauchy with Absorption	0.0061	0.5972	86.507
Cauchy, Absorp+Rough	0.0061	0.6921	86.506

The **Analysis Model** screen provides powerful features to analyze and visualize the Film Sense ellipsometric data.

Parameter	Value	Fit	Min	Max	St Incr
Rough(nm):	3.50	<input checked="" type="checkbox"/>	0	50	-1
Thick(nm).1:	87.12	<input checked="" type="checkbox"/>	5	500	10
n(633nm).1:	2.4491	<input checked="" type="checkbox"/>	0	6	-1
Bn.1:	0.010401	<input checked="" type="checkbox"/>	-0.1	0.1	-1
Cn.1:	0.0090553	<input checked="" type="checkbox"/>	-0.05	0.05	-1
Dn.1:	0	<input type="checkbox"/>	0	0	0
k(633nm).1:	0	<input type="checkbox"/>	0	0	0
k_slope.1:	0	<input type="checkbox"/>	0	0	0
%Grade.1:	7.39	<input checked="" type="checkbox"/>	-50	50	-1


## Features

- Standard, In Situ Multi-Layer, Multi-Sample, Trajectory, and Near Surface data analysis modes.
- Up to 10 model layers, with optional surface roughness, and substrate backside correction.
- Parameter ranges and starting increments to improve fit parameter convergence.
- Bruggeman effective medium approximation for mixed materials, and graded index layers.
- Cauchy, Sellmeier, Lorentz, Drude, Tauc-Lorentz, and Multi-Osc dispersion models.
- Temperature or composition dependent optical constant library files.
- Depolarization or transmission intensity data can be combined with the multi-wavelength data analysis.
- Simulate single measurement or dynamic data, and plot the Fit Diff vs. parameter value.
- Display fit parameter 90% confidence limits and correlation matrix, and estimate parameter accuracy


# Liquid Cell Option

## Powerful New Ellipsometric Characterization Capabilities


- Ellipsometric measurements in air and liquid ambients (“immersion ellipsometry”) are combined to provide more information and enhanced sensitivity for thin film characterization:
  - **Accurate index of refraction and thickness for very thin films:**  
≈ 0.01 index accuracy for native oxide film < 2 nm !
  - **n & k and Thickness for thin films, without a dispersion model**
- Easy to use and fill, preferred liquid is water (non-toxic, readily available)
- Film Sense software supports the analysis of the combined data sets
- Compatible with most Film Sense Ellipsometer systems (ex situ only, may not be suitable for all samples)



## Ellipsometric Measurements in Air and Water Ambients



## Combined Analysis: Air (CA) + Liquid (CL) Data



## ABOUT US

Our mission is to create compact and affordable multi-wavelength ellipsometer systems that provide optimal measurement solutions for a wide range of thin film applications.

Film Sense is located in Lincoln, Nebraska and was founded in 2013 by Mr. Blaine Johs, who has more than 30 years of experience in the ellipsometry industry.

Visit [www.film-sense.com](http://www.film-sense.com) for more information:  
*Customer Testimonials, Videos, Publication Lists, and more*

email us at [sales@film-sense.com](mailto:sales@film-sense.com)

