



L. Fink<sup>(1)</sup>, D. Podgorski<sup>(1)</sup>, M.U. Schmidt<sup>(1)</sup> and T. Hartmann<sup>(2)</sup>

A Stoe Stadi P powder diffractometer with Ge(111) monochromator yielding pure  $K_{\alpha 1}$ -radiation and the Dectris MYTHEN 1K detector has been chosen for PDF calculation experiments on Naphtalen ( $C_{10}H_8$ ).

For Cu  $K_{\alpha 1}$ -radiation the Stoe Stadi P has been equipped with a Dectris MYTHEN 1K with 320  $\mu$ m, for Mo  $K_{\alpha 1}$ -radiation with a MYTHEN 1K with 450  $\mu$ m and for Ag  $K_{\alpha 1}$ -radiation with a MYTHEN 1K with 1mm chip size. Synchrotron data has been taken at beamline X17A at NSLS Brookhaven ( $\lambda$ =0.1839Å). PDF calculations calculated with PDFgetX3 [1] yield a  $Q_{(obs\ max)}$  of 7.0 Å<sup>-1</sup> for Cu-, 11.4 Å<sup>-1</sup> for Mo-, 13.1 Å<sup>-1</sup> for Ag-K  $_{\alpha 1}$ -radiation and 19.5 Å<sup>-1</sup> for the synchrotron.



The direct comparison of the PDF curves of the synchrotron (yellow) and the Ag- $K_{\alpha 1}$ - experiment (blue), shows that the resolution of the Ag-data is amazingly similar!



Comparison of the laboratory setup (blue) and the synchrotron experiment (yellow)

Taking into account that  $\lambda_{(synchrotron)}$  has been app. ½ of Ag-K  $_{\alpha 1}$ , the measuring time (Ag-experiment 18h, synchrotron ½h) is more than reasonable for a laboratory setup.

This makes the Stoe Stadi P with Ag-tube and Dectris MYTHEN 1K an impressive alternative to expensive synchrotron experiments.

[1] Juhas, P., Davis, T., Farrow, C.L. and Billinge, S.J.L., J. Appl. Cryst. (2013). 46, 560-566.

<sup>(1)</sup> Institute of Inorganic and Analytical Chemistry, Goethe-University Frankfurt, Germany.

<sup>(2)</sup> Stoe & Cie GmbH Darmstadt, Germany