QUALITATIVE PDF ANALYSIS OF PIGMENT YELLOW 213¹⁾, $C_{23}H_{21}O_9N_5$ USING A STOE STADI P LABORATORY DIFFRACTOMETER WITH CU-, MO- OR AG- $K_{\alpha 1}$ RADIATION FROM SEALED TUBES AND A DECTRIS MYTHEN 1K DETECTOR

M. U. Schmidt, L. Fink, Institut f. Anorg. und Analyt. Chemie, Goethe Universität, Frankfurt/Main S. Correll, T. Hartmann, Stoe & Cie GmbH, Darmstadt


SETUP & RESULTS

To observe the influence of the wavelength and the detector noise using a laboratory powder diffractometer for the PDF determination and the maximum evaluable Q-value a 1mm glass capillary has been filled with Pigment Yellow 213 ¹⁾ (PY213) and adjusted on a goniometer head of a Stoe Stadi P powder diffractometer in Debye-Scherrer geometry at the Stoe & Cie application laboratory.


Picture 1: Pigment Yellow 213, C₂₃H₂₁O₉N₅

Data has been taken from set-ups with Cu-, Mo- and Ag-anode sealed tubes and the respective Ge(111)-monochromator for pure $K\alpha_1$ -radiation as well as a Dectris MYTHEN 1K detector with a wavelength optimised chip (320, 450 or 1000 μ m thicknes), a linear and an IP-PSD.

Data evaluation has been carried out by the Dept. of Inorg. Chemistry, Goethe University Frankfurt, using PDFgetX2 ²⁾.




Picture 2: Background and powder pattern of PY213 measured with Cu- and lin. PSD, Moand IP PSD as well as Agradiation and MYTHEN 1K (large to small)



Picture 3: F(Q) for PY213 with Ag-radiation

A Fourier transformation of the diffraction data F(Q) with $Q=4\partial(\sin\theta/\lambda)$ yields the observed PDF G(r), which shows the probability to find a pair of atoms in the distance r. The maximum theoretical Q value can be directly calculated from 2θ (max) of the goniometer.

For the Stadi P Q(theo max) could be 20 Å⁻¹ for Ag-radiation, but the observed Q(max) is much lower. G(r) yields a Q(obs max) for Cu-radiation of appr. 7 Å⁻¹, 10 Å⁻¹ for Mo- and 13 Å⁻¹ for Agradiation. Q(obs max) for ESRF Synchrotron data with λ = 0.40 Å has been 14 Å⁻¹3).



Picture 4: G(r) for PY213 measured with Cu-, Mo- and Ag-radiation (left to right)



Picture 5: Stoe Stadi P with MYTHEN 1K

As expected, the shortest wavelength yields the highest Q (obs max).

Because of the outstanding resolution and the low noise the Dectris MYTHEN 1K is the detector of choice and with the pure Ag $K_{\alpha 1}$ -radiation obtained by the Ge(111) monochromator the Stoe Stadi P powder diffractometer scores as well in XRD as in PDF calculation.

¹⁾ Schmidt, M.U. et al., Acta Cryst. (2009). B65, 189-199

²⁾ Qiu, X., Thompson, J.W. & Billinge, S.J.L. J. Appl. Cryst. (2004). 37, 678.

³⁾ Presentation of Martin U. Schmidt, EPDIC, Warszawa, Sept. 2008